

# Blood products, medicinal products or cultures need specific temperature storage

Newsteo data loggers enable to constantly control the temperature of climatic chambers and be alerted in real-time if there is any problem.



#### **ONLINE MONITORING**

WebMonitor web application enables any allowed user to bring up on screen measures from Newsteo loggers from any device connected on the internet.

## LOCAL MONITORING

RF Monitor Software enables to remotely configure and monitor all Newsteo products from a local PC.

## **©** CREATE YOUR OWN WIRELESS DATA LOGGERS NETWORK

Newsteo provides a wide range of stand-alone data loggers for any kind of fridge, freezer, climatic chamber....



#### **TEMPERATURE & HUMIDITY**

LOG22 [-40°C ... +70°C] [0 ... 100 %HR]


LGR64 [-40°C ... +85°C] [0 ... 100 %HR] IP65 Protected





## TEMPERATURE

LOM16 [-40°C ... +70°C] [0 ... 100%HR] LOG36 [-40°C ... +85°C] IP65 Protected EN12830




# ⊘ INSIDE PARCELS

USB mini loggers with automatic PDF report generation can be placed inside parcels to discreetly control the cold chain.

#### KIT FOR CHARACTERIZATION OF CLIMATIC AND THERMOSTATIC CHAMBERS FD X 15-140

This system allows you to characterize and qualify your enclosures such as climatic chambers, refrigerators, freezers. Depending on the size of the chamber, a 9 points or 15 points verification may be proposed.



# OUR CLIENTS

